



Service Manual

Lifecard CF

and

Lifecard 12

Firmware revision 7



SPACELABS
Healthcare

This Page is intentionally Blank

Contents Page

1.	PRODUCT OVERVIEW	7
1.1	INTRODUCTION	7
1.2	Specification	9
1.2.1	Standard Recording Mode – ECG Inputs	9
1.2.2	Standard Recording Mode - Pacemaker Pulse Detection	10
1.2.3	Standard Recording Mode - Data Storage	10
1.2.4	User Interface	10
1.2.5	Power Requirements	11
1.2.6	Extended Recording Mode	11
1.2.7	Additional Specifications in 12-Lead Mode	12
1.2.8	Physical and Environmental	12
1.2.9	Electro-Magnetic Compatibility	13
2.	SAFETY AND REGULATORY	15
2.1	Intended Use of Equipment	15
2.2	Safety Classification	15
2.5	Explanation of Markings	17
3.	LCF BUILT-IN TESTS	21
3.1	Primary Built-in Tests	21
3.2	Second Level Built-in Tests	21
3.3	Third Level Built-in Tests	21
3.4	User Level Tests	22
3.5	Fault Protection While Recording	22
3.6	Lifecard CF Menu Options	23
3.6.1	Main Menu Options	23
3.6.2	Setup Menu	25
3.7	Lifecard CF Sounds	28
3.8	Error and Warnings Displays	29
3.9	Decontamination	32
3.10	Patient Cable	33
3.11	Battery	34
3.12	Battery Check	34
3.13	Checking the Hardware & Software Revision	34
3.14	Compact Flashcard	35

4.	Test Procedure	37
4.1	Checking the Patient Cable	37
4.2	LC12 Varios Yoke Assembly	37
4.3	Lifecard CF Check-out	39
5.	Spares & Part Numbers	41
6.	Assembly Views	43
6.1	Front Assembly - Front View - Complete	43
6.2	Front Assembly - Rear View	44
6.3	Intermediate Moulding	45
6.4	Circuit Board	46
6.5	Patient Cable Mounting with Clip	48
6.6	Short Patient Cable	49
6.7	Varios with 46-1123 Cable 46-1152	50
6.8	Varios Lead Parts	51

©*Copyright 2010 Spacelabs Healthcare Ltd.*
1 Harforde Court, John Tate Road, Hertford. SG13 7NW

This Page is intentionally Blank

1. PRODUCT OVERVIEW

1.1 Introduction

The Lifecard CF is a compact Holter Ambulatory ECG Recorder utilising a digital storage technique to store the ECG recording onto a Compact Flash (CF) card. The Lifecard CF provides continuous recording of 2 or 3 leads of ECG for up to 48 hours in standard mode and up to 7 days in extended mode.

The Lifecard 12 option provides continuous recording of 12 leads of ECG for a period of 24 hours. The recorder has a built in display for you to monitor the ECG and pacing detection during hook-up. This enables you to verify the ECG quality before starting the recording. Menu options are selected using the 2 buttons on the front of the recorder unit.

The Lifecard CF requires one AAA battery. The patient cables for the Lifecard CF are designed to prevent accidental disconnection from the recorder by the patient.

The Patient Event button on the front of the recorder unit enables the patient to indicate symptomatic episodes in the recording for correlation with the patient diary. Pacemaker pulse detection may be enabled and disabled by the physician or cardiac technician. Recordings may be analysed using a Spacelabs Healthcare Pathfinder, Impresario, or Lifescreen Holter analysis system, if they have compatible hardware and software. (Lifescreen is incompatible with 12-lead recordings.)


The Lifecard CF comprises two sections, the 'Recorder Unit' and the 'Patient Cable Unit'.


LIFECARD CF RECORDER UNIT

Recorder Unit



Patient Cable Unit





1.2 Specification

1.2.1 Standard Recording Mode – ECG Inputs

Channels	3, type BF applied part patient isolation types 2 or 3 channel 3 electrode, 2 4 electrode, 3 channel 6 electrode with detachable leadwires.
Input impedance	> 5M Ω Ohms
Input	DC offset \pm 300 mV, with saturation recovery circuit (3 seconds max)
CMRR	> 60 dB at 10Hz, > 80 dB at 50 Hz and above, 2 Vpp signal
Dynamic range	10 mV
Resolution	2.5 μ V

Calibration	± 5%
Bandwidth	0.05 - 40 Hz (-3 dB)
Sampling rate	1024 samples per second per channel
Noise filter	Linear phase filter effective from 60Hz to > 1MHz, 128 samples per second output rate

1.2.2 Standard Recording Mode - Pacemaker Pulse Detection

Sensitivity	7 mV nominal, channels 1 and 2 only
Noise rejection	> 50 mVpp for sinusoids up to 200 Hz
CMRR	2 V common mode spikes are rejected
Refractory period	40 ms

1.2.3 Standard Recording Mode - Data Storage

Media type	Removable card, CompactFlash Association standard (Type 1)
Data types	Full disclosure ECG, with pacing and patient event markers. Recording Time and Date. Patient name and record number (Pathfinder systems), Encrypted patient record file (CardioNavigator). 8 second voice recording. Recorder serial number
Capacity req	15 Mbytes per channel per 24 hours eg. a 48 hour three channel recording occupies 90 Mbytes

1.2.4 User Interface

Type	Text menus with audio cues and keys for up, down and select
Languages	English, German, French, Italian, Spanish, Danish and Polish languages also Hebrew patient ID support
Clock	Clock and calendar (to 2098), selectable 12/24 hour and US/European date formats. 13mm digit height for patient use

Basic features	Pacing detection on/off, hook-up display, voice recording for patient identification
Ancillary features	Identify and delete unread recordings, warning/error screens for battery and memory card conditions
Hook-up display	Real time display of each channel, with 60 μ V/30 ms resolution and pacing annotation
Set-up options	Time and date, language, display contrast, recorder identification

1.2.5 Power Requirements

Disposable cell	Single AAA alkaline (Duracell MN2400 or equivalent), two 24 hour recordings or one 48 hour recording
Rechargeable cell	Single AAA nickel metal hydride (Ansmann 600 mAh or equivalent), one 24 hour recording per charge
Battery check	User is warned of poor battery condition before recording
Clock battery	Internal rechargeable cell, charged during recording. The clock is maintained for > 3 months between recordings

1.2.6 Extended Recording Mode

Channels	2 channel recording, with pacing detection
Cable Types	2 channel 3 electrode or 2 channel 4 electrode
Resolution	10 μ V
Sampling Rate	256 samples per second per channel
Compression	10 μ V maximum compression error when tested with MIT-BIH Arrhythmia and Compression databases
Capacity required	90 Mbytes for dual channel 1 week
Disposable Cell	AAA alkaline (Duracell MN2400 or equivalent) for 1 week.
User Interface	Includes an audible alarm to alert the patient if an electrode becomes detached. Sense current is < 10 nA.

Note: other specifications are the same as Standard Recording Mode.

1.2.7 Additional Specifications in 12-Lead Mode

Channels	Standard 12-lead, one neutral and nine active electrodes
Cable types	10 electrode, defibrillation protected, IEC or AHA code
Isolation	DEFIBRILLATION-PROOF TYPE CF APPLIED PART
Input impedance	10 Mohm
Sampling rate	4096 samples per second per channel
CMRR	> 80dB per IEC and ANSI/AAMI methods
Suppression	Active neutral system ('right leg drive')
Resolution	0.6 μ V
Noise	< 0.6 μ V RMS
Pacing detection	>2mV, 200 μ s to 5ms pulse in any electrode
Capacity required	256 MByte card for 24 hour recording
Battery	An alkaline AAA cell is required for 24 hour recording
Fault tolerance	In the event of electrode detachment noise is suppressed, and the available leads are recorded (differential V leads only if R, L or F is detached)

Note: other specifications are the same as Standard Recording Mode

1.2.8 Physical and Environmental

Dimensions	96 x 57 x 17.5mm with patient cable fitted
Weight	Recorder body 55g: patient load 130g including battery, card and typical patient cable
User labelling	Area provided is 52 x 15 mm
Temperature	0 to 45°C operation, -20 to 65°C storage
Humidity	Operation or storage 5% to 95%, non-condensing
Pressure	Operation or storage air pressure 700 - 1060 mbar
Shock	1 m drop

1.2.9 Electro-Magnetic Compatibility

General	<i>Complies with EN60601-1-2:1993 and ANSI/AAMI EC38:1998</i>
ESD (1)	<i>4 kV air and 2 kV contact discharges: no interruption in recorder function</i>
ESD (2)	<i>8 kV air and 6 kV contact discharges: no damage to the recorder, recording resumes automatically in < 10 s</i>
Radiated emissions	<i>CISPR 11:1997, EN55011:1998 Group 1 Class B</i>
Radiated immunity	<i>3 V/m 26 MHz - 1 GHz, 80% AM modulated at 5 Hz. Keyed carrier immunity to EN50082:1996</i>

This equipment has been tested and found to comply with the limits for a class B computing device in accordance with the specifications in Subpart J of Part 15 of FCC Rules, which are designed to provide reasonable protection against interference to radio and television reception.

This equipment generates and uses radio frequency energy and if not installed and used in accordance with the instructions it may cause interference. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause interference to radio or television reception, which can be determined by turning the equipment off or on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient the receiving antenna
- Relocate the equipment with respect to the receiver
- Move the equipment away from the receiver

If necessary, the user should consult Spacelabs Healthcare or an experienced radio/television technician for additional suggestions. The user may find the following booklet prepared by the Federal Communications Commission helpful:

“How to Identify and Resolve Radio-TV Interference Problems”

This booklet is available from the U.S. Government Printing Office, Washington, DC 20402, Stock No. 004-000-00345-4.

This Page is intentionally Blank

2. SAFETY AND REGULATORY

2.1 Intended Use of Equipment

The Lifecard CF Holter recorder is to be used for the non-invasive ambulatory recording of two or three channel electrocardiograms on a standard commercial compact flash card.

The Lifecard 12 option is to be used for the non-invasive ambulatory recording of 12-lead electrocardiograms on a standard commercial compact flash card.

The recorder allows data to be collected over a continuous period of up to 7 days whilst allowing the subject to perform most of their normal daily activities.

The recordings can be analysed on compatible analysis systems from Spacelabs Healthcare.

This device has been designed and supplied specifically for the long term recording of electrocardiograms in ambulatory patients using standard Holter monitoring techniques. It shall not be used for any other purposes.

The device shall be operated only by suitably competent personnel trained in the use and procedures of Holter electrocardiography for diagnostic purposes.

The Lifecard CF comprises two sections; the 'Recorder Unit' and the 'Patient Cable Unit'.

2.2 Safety Classification

This device has been designed in accordance with EN60601 - 1 , "Medical electrical equipment, Part 1: General requirements for safety", as follows:

1. EQUIPMENT with an INTERNAL ELECTRICAL POWER SOURCE. The equipment is designed to be battery operated only. Under NO circumstances shall a mains powered battery eliminator or any other external power source be used with the equipment.

-
2. EQUIPMENT having a TYPE BF APPLIED PART. or
 3. EQUIPMENT having a TYPE CF APPLIED PART if so marked.
 4. IPX4 EQUIPMENT protected against the ingress of splashing water, if so marked. Otherwise, ORDINARY EQUIPMENT, without protection against ingress of liquid.
 5. Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide, or flammable cleaning agents.
 6. Rated for CONTINUOUS OPERATION.
 7. EQUIPMENT with an APPLIED PART, specifically designed for application where a CONDUCTIVE CONNECTION is made to the PATIENT, but not directly to the heart. According to ANSI/AAMI EC38:1998.Lifecard CF is Type 1 ambulatory ECG device.

2.3 Adjustment, replacement of parts, maintenance and repair

The device requires no routine adjustments to maintain its operation.

The device contains no user serviceable parts. It shall be serviced only by Spacelabs Healthcare or by an agent accredited by them to service device of this type. Unauthorised repairs or dismantling of the device will invalidate the guarantee.

2.4 Defects and abnormal stresses








For continued safety the device must not be maltreated, used outside its specified operation conditions, or stored outside its specified storage conditions.



Lifecard CF contains protection against electrostatic discharge, but there is no protection against defibrillators. To avoid damage the device should be removed before defibrillating. The Varios active yoke and 46-1123 / 46-1127 patient cables have defibrillator protection. (The protection is a combination of the cable and the yoke).

Whenever it is likely that protection has been impaired, the device shall be made inoperative and secured against any unintended operation. The protection is likely to be impaired if, for example, the device shows visible damage.

- a) shows visible damage
- b) fails to perform the intended measurements
- c) has been subjected to prolonged storage under unfavorable conditions
- d) has been subjected to severe transport stresses
- e) the device has been connected to a patient during defibrillation.

2.5 Explanation of Markings

	CE Mark
	Date of Manufacture
	Manufacturer
	Consult Documents This symbol means you should read the accompanying documents
	EQUIPMENT having a TYPE BF APPLIED PART
	EQUIPMENT having a DEFIBRILLATION-PROOF TYPE CF APPLIED PART
	Battery Eject symbol

	Protected against ingress of splashing water when the patient cable is fitted.
AECG-	Type 1 ambulatory ECG device according to ANSI/AAMI EC38:1998
	This product must be sent to separate collection facilities for recovery and recycling
12 LEAD ✓	Compatible with Varios Active Yoke & 10-electrode cable

2.6 Warranty

Subject to the conditions set out below, Spacelabs Healthcare (“The Company”) warrants that its Products will be free from defects in material and workmanship for a period of 12 months from delivery.

This warranty is given by The Company subject to the following conditions:

1. The Company shall be under no liability in respect of any defect arising from fair wear and tear, willful damage, negligence, abnormal working conditions, failure to follow instructions (whether oral or in writing), misuse, improper installation or alteration or repair of the Products without The Company’s approval.
2. The above warranty does not extend to parts, materials or devices not manufactured by The Company, in respect of which the Customer shall only be entitled to the benefit of any such warranty or guarantee as is given by the manufacturer to The Company.
3. Subject as expressly provided here, all warranties, conditions or other terms implied by statute or common law are excluded to the fullest extent permitted by law.
4. Any claim by the Customer which is based on any defect in material or workmanship of the Products shall be notified to The Company immediately after discovery of the defect. If the Customer does not notify The Company accordingly, the Customer shall not be entitled to reject the Products and The Company shall have no liability for such defect.
5. Where any valid claim in respect of any of the Products which is based on any defect in the material or workmanship of the Products is notified to The Company, The Company shall be entitled to replace or repair (at The Company’s sole discretion, either at the Customer’s premises or at The Company’s premises in the United Kingdom) the Products (or part in question) but The Company shall have no further liability to the Customer.

-
6. The Company shall not be liable to the Customer by reason of any representation, or implied warranty, condition or other term, or any duty at common law, or for any consequential loss or damage (whether for loss of profit or otherwise), costs, expenses or other claims for consequential compensation whatsoever arising out of or in connection with any act or omission of The Company relating to the manufacture or supply of the Products or use by the Customer.
 7. Spacelabs Healthcare recommends the use only of approved accessories and parts. The use of third party accessories may result in damage to recordings or equipment, and may invalidate your warranty.

2.7 Contact Details



Spacelabs Healthcare Ltd
1 Harforde Court
John Tate Road
Hertford
SG13 7NW
United Kingdom

Tel: +44 (0)1992 507700
Fax: +44 (0)1992 501213

Spacelabs Healthcare Inc.
5150 220th Ave. SE
Issaquah,
WA 98029
USA


Tel: +1 425 657 7200
Fax: +1 425 657 7212

Web site: www.spacelabshealthcare.com


2.7 Modifications

For continued safety, the device shall not be subjected to any Unauthorised modifications and must be used only for the purpose for which it was originally supplied.

CIRCUIT BLOCK DIAGRAM



LIFECARD 12 BLOCK DIAGRAM



3. LCF BUILT-IN TESTS

3.1 Primary Built-in Tests

When a battery is inserted the internal functions of the microcontroller are tested automatically. If any errors are found the recorder will start to beep loudly, before the patient cable is fitted. The pattern of beeps indicates the cause of the problem.

3.2 Second Level Built-in Tests

When the patient cable is fitted the microcontroller applies further automatic tests before proceeding with normal operation. The integrity of the complete operating programme is tested at this point. If any errors are found the recorder will start to beep after the cable has been fitted. The pattern of beeps again indicates the cause of the problem.

3.3 Third Level Built-in Tests

Lastly, the microcontroller tests the remaining peripheral circuits, communicating any errors with a text message on the LCD. The following service screen codes are used:

- 0001 Oscillator/clock fault
- 0002 SRAM fault
- 0003 EEPROM data error (data includes the unit serial number, etc.)
- 0004 ADC operation fault
- 0005 ADC calibration fault
- 0006 Configuration fault
- 0007 LCD Voltage fault
- 0008 Varios Active Yoke EEPROM data error
- 012 (TIMI recorders only) Clock, EEPROM, SRAM or LCD Voltage fault.

3.4 User Level Tests

Before commencing the recording the recorder checks all aspects that are under the user's control, and will display warning or error messages as described in the user manual. Items checked include:

- Battery voltage
- Loss of time and date since last use
- CF card needs conditioning
- CF card contents- have they been read out for analysis?
- Patient cable type

If a recorder frequently requires the time and date to be entered or corrected, it should be returned for service.

3.5 Fault Protection While Recording

During recording the recorder operation is supervised in various ways, to ensure that all recorded data is correct and the recording is as complete as possible:

- The battery level is continually monitored. Low battery shutdown is described in the user manual.
- If the LCD is found to be in the wrong state the recorder will re-load the display.
- Data written to the CF card includes error detection information which is checked before analysis.
- The recorder shuts down immediately if the cable is removed.
- If the patient cable is removed and then replaced, the recorder will continue after a countdown message. The analyser will display a corresponding blank period in the ECG.
- The 32768Hz crystal is monitored via a watchdog to ensure that it is operating correctly.

Other transient faults can cause a short blank in the ECG while the recorder resets and restarts:

- Intermittent power supply, due to contact corrosion, bad joints, or possibly severe shock.
- The recorder software running improperly (watchdog timer).
- Error messages from the CF card that cannot be resolved.
- Any other unexpected condition such as data errors due to large electrostatic discharge.

3.6 Lifecard CF Menu Options

The Lifecard displays are controlled by two buttons:

- The yellow ▲ Up and ▼ Down button is used to move up and down the menus, to highlight the menu option. You also use it to choose which channel is displayed on the monitor.
- The green ► Select button is used to select the highlighted option.


3.6.1 Main Menu Options

Displays options for starting and setting up the recorder.




1. Lifecard CF Option

If alternative operating modes are available on your recorder and you wish to use them, contact Spacelabs Healthcare Ltd. for the appropriate user manual.




2. Pacing option

Use this option to select or deselect pacing detection.



3. Start... Option

Select Start when you want to make a recording in Standard Recording Mode, up to 48 hours of 2 or 3 channel ECG, or 24 hours of 12-lead ECG.



4. Start Week ...Option

Select Start Week when you want to make a recording in Extended Recording Mode, to record up to 7 days of continuous 2 channel ECG.



5. Recorder ID

This is a name or number you have chosen for the recorder using an option in the Set Up menu. Your chosen ID will appear in the Main Menu and on Pathfinder reports.

There is a space inside the recorder unit for you to affix a bar code or similar.

6. Set Up Option

Select this option to configure the recorder.

7. Language Option

Select this to change the language displayed on the monitor. Keep pressing the green ► Select button to switch between the languages until the correct one is displayed.


8. About option

This option displays the Software Version number and Hardware Serial number. Pressing the green ► Select button returns you to the Main Menu.


3.6.2 Setup Menu

1. Contrast Option

Displays a bar indicating the current setting of the contrast.




The display can be darkened or lightened by pressing the yellow Up or Down button. Press the green Select button to return to the Set Up Menu.



2. Time Adjustment Option

This option enables you to select a 12 or 24 hour clock display. Use the yellow ▲ Up or ▼ Down button to change the selection and then press the green ▶ Select button to select it.



Adjust the hours by pressing the yellow ▲ Up or ▼ Down button until the correct one is displayed, then press the green ► Select button. This moves the highlight onto the minutes which you then adjust in the same way.




Pressing the green ► Select button again returns you to the Set Up Menu.

3. Date Adjustment Option

This option enables you to select a European or American date format. Use the yellow ▲ Up or ▼ Down button to change the selection and then press the green ► Select button to select it.




Adjust the date in the same way as you have adjusted the time.




4. Speech Recording On/Off

Press the green ► Select button to switch between the Speech Recording option being on or off.



5. ID:Enabled/Disabled

This option enables you to give the recorder a personalised name or number ID. Your chosen ID will appear in the Main Menu and on reports.




To use this option, highlight the ID:Disabled line in the Set Up menu and press on the green ▶ Select button.

Press either of the yellow ▲ ▼ buttons until 'ID: enabled' is displayed.

Then press the green ▶ Select button to enter your recorder ID.


The first character of the ID will be highlighted. Use the yellow ▲ ▼ buttons to change the character as required and then press the green ▶ Select to move the highlight onto the next character.



You can enter a name or number up to 10 characters long.


6. Sounds On/Off

Use the green ▶ Select button to select or deselect the Lifecard sounds.



7. Service

The recorder will remind you when the annual service is due.



3.7 Lifecard CF Sounds

The Lifecard CF generates sounds to confirm your actions, or to inform you about the status of the recorder.

1. Click

The 'click' sound is emitted whenever you press the yellow ▲ Up or ▼ Down button.

2. OK 'boing'



The 'OK' boing tells you that the recorder is functioning properly and also sounds whenever you press the green ► Select button. The tone is a higher note than the 'error' tone.

3. Error tone



The Error tone tells you that the recorder has identified an error or warning condition. The tone is a lower note than the 'OK' boing.

4. High Impedance or Lead Off Alarm




During an extended recording, if a lead becomes detached, the Lifecard will generate an alarm to warn the patient. The alarm is a 'ding-dong' tone lasting 15 seconds, and the clock display will be replaced by the LEAD OFF message, indicating which channel has been affected.

Once the patient has replaced the lead, the clock display resumes.

You can turn the sound off by selecting Sound Off in the set up menu. This does not affect the 'Lead Off' alarm.

3.8 Error and Warnings Displays


1. You will see this error message if you are trying to use Extended Recording Mode with a six electrode cable. Either select Start. . . for Standard Recording Mode or change the patient cable.




2. You are trying to use Standard Recording Mode with a two electrode cable. Either select Start Week for Extended Recording Mode or change the patient cable.



3. There may be a fault in the cable, or you may require a firmware upgrade to use this cable type.




4. If this error message is displayed, contact the Service Department at Spacelabs Healthcare Ltd.




***ERROR 0007
Service required
by
Del Mar Reynolds

5. This message is displayed if a patient loses a lead during an Extended Recording. The Lifecard will also generate an alarm which will consist of a series of 'dingdong' tones lasting 15 seconds. Once the patient has replaced the lead, the clock display resumes.



***LEAD OFF ***
Channels 1 and 2
11.11.03 15:49
Recording


6. This warning is generated if you have selected Start Week... and the battery is partially discharged.



WARNING
Battery Low
1.36 V
▶ Continue

7 day recording needs a new alkaline battery.

7. You will see a warning of this format if you have selected Start Week... and the capacity is too small to record 7 days. Select Continue only if the recording description meets or exceeds your requirements.




WARNING
2 Ch. for 3 days
Capacity 45.8 MB
▶ Continue


8. This message appears if the card capacity and patient cable are not appropriate to provide a 24 hour recording.




9. The Lifecard has been closed without a Flashcard inserted. Reopen the recorder and insert the card.




10. The Flashcard has been inserted before being correctly initialised. Remove the card and initialise the card.




11. If the Flashcard contains a recording that has not previously been analysed, you have the option to delete it. Select the 'DETAILS' option - if the patient information is available, then the Patient screen is displayed. Otherwise the screen will display just the time that the recording was started.



12. If the recorder detects that the rechargeable clock battery has discharged since the last use, you must enter the time and date.



13. If this message is displayed before you start a recording, you should insert a new battery immediately.



3.9 Decontamination

Before commencing any service or maintenance procedures, ensure the Lifecard CF recorder has been suitably decontaminated. Unless contamination with body fluids is known or suspected, we recommend low-level disinfection.

See the user instructions for further details on cleaning and disinfecting the recorder.

In the case of severe contamination, the recorder may be beyond repair. Please contact Spacelabs Healthcare for further advice.

3.10 Patient Cable

Lifecard CF patient cables use screened lead-wires with tinsel conductors and an antimicrophonic barrier to provide strength, flexibility and low noise. There is a choice of 7 patient cable units:

- 4 electrode cable with belt clip
- 4 electrode cable, short, no belt clip
- 3 electrode cable with belt clip
- 3 electrode cable, short, no belt clip
- 6 electrode cable with belt clip and renewable lead wires (unscreened)
- 10 electrode cable for Varios Active Yoke, AHA code
- 10 electrode cable for Varios Active Yoke, IEC code

Two electrode cables for single channel ECG are no longer offered. We recommend three electrode hook-up in long term applications.

Cables with 4 electrodes provide two channels of ECG. The popper colour code, which is detailed on the cable, conforms to AHA (1985) recommendations. The short cable without belt clip is designed to maximise patient comfort when the recorder is worn under clothes.

The 3 electrode scheme and popper colour code, as detailed on the cable label, is proprietary to Spacelabs Healthcare.

The 6 electrode cable provides three channels of ECG from 6 electrodes. The individual lead wires are detachable and can be replaced with any DIN42802 connector lead wire,

useful in situations requiring unusual hook-up or electrode terminations. The yoke is colour coded to AHA (1985) recommendations.

The cable yoke/back unit forms a complete assembly and is not serviceable. In the case of failure, the complete assembly must be replaced.

In the case of the 6 electrode cable, the individual lead wires may be replaced if broken.

The 10 electrode cable is replaceable as a complete assembly by releasing the screw securing it to the Varios active yoke.

The Active Yoke itself is not serviceable and must be returned to Spacelabs Healthcare for repair.

3.11 Battery

Lifecard CF requires one AAA cell, either:

Disposable cell	Single AAA alkaline (Duracell MN2400 or equivalent), two 24 hour recordings or one 48 hour recording or in Extended Recording Mode up to 7 days, or in 12-lead mode 24 hours.
Rechargeable cell	Single AAA nickel metal hydride (Ansmann 600 mAh or equivalent), one 24 hour recording per charge, or in Extended Recording Mode up to 3 days. Not suitable for 12-lead recording

3.12 Battery Check

The user is warned of poor battery condition before recording, when starting the recorder. Please refer to the detailed instructions, in this user manual.

The real-time clock is maintained by an internal rechargeable lithium cell, charged during recording from the main battery. With a full charge, the clock is maintained for at least 3 months after the main battery is removed or exhausted.

The clock cell is not replaceable by the user, and in the case of suspected failure the Lifecard CF should be returned to Spacelabs Healthcare for service.

Dispose of used batteries carefully, using environmentally friendly methods wherever possible.

3.13 Checking the Hardware & Software Revision

1. Insert a battery and close the back cover.
2. Press the Yellow (DOWN) key to highlight the About menu option, then press the Green (SELECT) key.
3. The revision numbers are displayed on the LCD.

3.14 Compact Flashcard

The recorder uses a CF card (type 1) Compact Flashcard, non volatile memory card meeting the Compact Flash Association type 1 standard.

A new flashcard must be conditioned before first-time use, which may be done by the user from a menu option in the Spacelabs Healthcare Ltd. Holter analyser. However, cards supplied by Spacelabs Healthcare are already conditioned.

Full details are in the user instructions manual. If the card has not been supplied by Spacelabs Healthcare Ltd, ensure the type is exactly as listed in the Accessories section of this manual Other types may give short or unreadable recordings.

This Page is intentionally Blank

4. Test Procedure

4.1 Checking the Patient Cable

- Fit a battery and CF card to the Lifecard CF
- Connect the patient cable to the 9270 Tracker calibrator. (For a 3 lead cable – Red to CH1+, Yellow to CH2-, Green to CH2+, and select Cal1, PACE 2.) If using another cable or source ensure that there is a unique signal in each channel.
- In the main menu, select Start....
- Skip the voice recording
- Verify that the patients name has been written on the card.
- In the hook-up screen, scroll down the three channels and check signal and quality when the cable is gently flexed to reveal any intermittent cable faults. Channel 2 has a negative pacing spike only.
- When all the leads have been checked, disconnect the Lifecard CF and calibrator.

4.2 LC12 Varios Yoke Assembly


The Lifecard CF used to check these leads must have version V07.xx firmware. An Electromedicina ST-20 ECG Simulator is required, select Pacemakers/DVI, test signal.

- Fit the battery and CF card to the Lifecard CF
- Connect the patient cable to the ST-20 simulator:

ST-20	N	R	1	2	3	4	5	6	L	F
H10EU cable	N	R	C1	C2	C3	C4	C5	C6	L	F
H10US cable	RL	RA	V1	V2	V3	V4	V5	V6	LA	LL

- Connect the Lifecard CF to the Varios yoke
- Verify that the Lifecard CF shows the Lifecard 12 splashscreen
- In the main menu, select Pacing On
- In the main menu, select Start...
- Skip the voice recording
- Verify that the patients name has been written on the card.
- In the hook-up screen, scroll down through 9 leads, checking each one:
- All should have regular pacing markers (small vertical arrow with P).

-
- All should have baselines that settle in the central 1/3 of the screen.
 - The expected signals are shown in the figure below.
 - Checking leads I, II and III is not necessary.
 - When all the leads have been checked, disconnect the Lifecard CF and ST-20 simulator.



4.3 Lifecard CF Check-out

Connect the 3 Lead Patient Cable to the 9270 Tracker Calibrator, Red to CH1+, Yellow to CH2-, Green to CH2+, and select Cal 1, PACE 2. If using another cable or source ensure that there is a unique signal in each channel.

Fit a CF card, a AAA battery and the 3 Lead Patient Cable to the unit under test. Check that the speaker tone and display are normal. In the main menu ensure that PACING DET. Is ON, and then select START.

Select RECORD and speak clearly into the microphone from a distance of approximately 150 mm. The recording is replayed automatically; check that the recording / replay is clear and audible.

Select CONFIRM and then check the Channel 1, Channel 2 and Channel 3 hook-up displays (push the 'down' button twice). Using the 9270 calibrator, Channel 2 has a negative pacing spike only.

Press the 'select' button to leave the hook-up display then select START NOW.

Approximately 30 minutes after starting the recording press the Green button to mark a Patient

Event in the recording. Wait for 1 minute and remove the Patient Cable. Wait for another minute before re-attaching the Patient Cable. After 30 seconds the Recording should continue.

Leave the recording running for up to 24 hours.

Remove the 3 Lead Patient Cable and eject the CF card, eject the AAA battery to check for jamming.

Analyse the recording on Pathfinder or Lifescreen. Check for the Patient Event and recording break 30 minutes into the recording.

CHECK LIST

Hospital:

**Fault Reported:
Software Version before
updating:**

Serial No:

		Checked	Replaced
1	Front case damage	Cracks in moulding - Replace front assy. Lens scratched / fogging	
		<input type="checkbox"/>	<input type="checkbox"/>
2	Patient Cable	Cracks in moulding - Replace complete cable Catch moulding Hinge moulding Belt clip Connecting pins	
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
3	Battery contacts	Replace if damaged	
		<input type="checkbox"/>	<input type="checkbox"/>
4	Insert firmware update card + battery	Connect patient cable Unit will bleep for 10secs to update 1. company logo will appear 2. Lifecard 7 day will appear if not 7 day - select from Lifecard heading	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
5	Display	Check for missing lines	
		<input type="checkbox"/>	<input type="checkbox"/>
6	Select Setup menu	Contrast (adjust if necessary) Set Time Set Date Set speech Rec On Set sounds On Select "Back" to exit set up menu	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
7	Current software version	"About" heading note S/W version upgrade via CF card if necessary	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
9	Insert initialised flashcard (connect patient cable)	Enter "Start" make voice recording - confirm check for distortion	
		<input type="checkbox"/>	
10	Check leads for fractures (connect to calibrator)	Monitor signal on display Fracture at grommet Fracture at lead separation	
		<input type="checkbox"/>	
		<input type="checkbox"/>	
11	Test recording	Up to 24 Hours	
		<input type="checkbox"/>	
12	Test results on Pathfinder	Print results if required	
		<input type="checkbox"/>	


5. Spares & Part Numbers

Product Code	Description
19-5000	Lithium Battery, Rechargeable
19-5044	LCF Speaker, 64 Ohms
19-5074	Microphone, Condenser
19-5075	Battery Contact, Positive
19-5077	LCD Display
19-7508	Screw PT 1.8 x 10mm, Philips
23-0223	Pocket Clip
23-0224	Mount
23-0226	Spring Pocket Clip
46-0037	Hinge Moulding
46-0043	Top Cover Moulding
46-0044	Intermediate Moulding
46-0045	Card Ejector Moulding
46-0046	Retainer Moulding, Up/Down
46-0047	Retainer Moulding, Select
46-0050	Switch Moulding, Up/Down
46-0051	Switch Moulding, Select
46-0056	PCB Insulator
46-0057	Internal Label
46-0059	Main Seal
46-0062	Top Cover Assembly
46-1570	LCF PCB Assembly, Exchange
46-0420	Mic Speaker
46-0421	Speaker Seal
46-0426	Acoustic Membrane
46-0462	Battery Contact, Negative
46-0493	Catch Moulding Printed
46-0551	Waterproof Lens
46-0418	3 Electrode Patient Cable with Clip
46-0480	3 Electrode Patient Cable, Short


46-0556	4 Electrode Patient Cable with Clip
46-0557	4 Electrode Patient Cable, Short
046-0031	6 Electrode Patient Cable with Clip
19-7509	6 Electrode Wire Lead Set
046-0523	7 Electrode Patient Cable with Clip
19-7510	7 Electrode Wire Lead Set
46-1152	LC12 Varios with 46-1123 Cable
46-1153	LC12 Varios with 46-1127 Cable
46-1122	LC12 US Patient Cable (Leads)
46-1123	LC12 EU Patient Cable (Leads)
46-1125	LC12 Varios Active Yoke (No Leads)

6. Assembly Views


6.1 Front Assembly - Front View - Complete




6.2 Front Assembly - Rear View




6.3 Intermediate Moulding




6.4 Circuit Board







6.5 Patient Cable Mounting with Clip



6.6 Short Patient Cable




6.7 Varios with 46-1123 Cable 46-1152




6.8 Varios Lead Parts

LC12 EU PATIENT CABLE 46-1123



VARIOS ACTIVE YOKE (NO LEADS) 46-1125



www.spacelabshealthcare.com

Spacelabs Healthcare Inc.
5150 220th Ave. SE
Issaquah
WA 98029, USA

Tel: +1 425 657 7200
Fax: +1 425 657 7212



Spacelabs Healthcare Ltd.
1 Harforde Court
John Tate Road
Hertford, SG13 7NW, UK

Tel: +44 (0) 1992 507700
Fax: +44 (0) 1992 501213

R_x only



CE
0120